# TriSensor

For product code: KIT035 - Incubation temperature: 40°C

GB

**A Rapid Test For  $\beta$ -lactams,  
Sulfonamides and Tetracyclines Detection in MILK**

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
**Figure A : Visual interpretation**

**TABLE A : Limits of detection of Trisensor**

B-LACTAMES (ppb)		SULFONAMIDES (ppb)		TETRACYCLINES (ppb)	
<b>Penicillines</b>		Sulfadiazine	8-10	Tétracycline	80-100
Penicilline G	2,5-3,5	Sulfapyridine	0,5-1	Oxytétracycline	60-70
Ampicillin	3-4	Sulfathiazole	7,5-8,5	Chlortétracycline	50-60
Amoxicillin	3-4	Sulfamethoxazole	320-360	Doxycycline	20-30
Oxacillin	13-15	Sulfamethazine	1-2		
Cloxacillin	7-9	Sulfamethoxyipyridazine	2-3		
Dicloxacillin	4-6	Sulfadimethoxine (Sdm)	10-15		
Nafcillin	50-70	Sulfacetamide	300-600		
<b>Cefalosporines</b>		Sulfamerazine	2-3		
Ceftiofur	8-12	Sulfamonomethoxine	8-12		
Cefquinome	16-20	Sulfaquinoxaline	14-18		
Cefazolin	15-19	Sulfachloropyridazine	5-10		
Cefapirine	4-6	Sulfaguanidine	15-25		
Cefacetile	16-20	Sulfamethizole	220-260		
Cefoperazone	3-4	Sulfasalazine	250-350		
Cefalexin	700-800				
Cefalonium	3-5				

Trisensor is a rapid test that allows you to simultaneously detect the presence of both  $\beta$ -lactam, Sulfonamide and Tetracycline molecules in a milk sample.

## II. Summary of the protocol

- Add 200  $\mu$ l of milk into one reagent microwell and mix to homogeneity;
- Incubate 3 min at 40°C;
- Dip one Dipstick into each microwell;
- Continue incubating for 3 min at 40°C;
- Read the colour intensities.

## III. Reaction Mechanism

Trisensor is a competitive test involving two receptors and generic monoclonal antibodies in one single operation. The test requires the use of two components. The first component is a microwell containing predetermined amounts of both receptors and antibodies linked to gold particles. The second is a dipstick made up of a set of membranes with specific capture lines. For a valid test, the upper red control line has to be visible after the second incubation (see Figure A). The other three are the specific "test" lines placed below the control line. The line for  $\beta$ -lactam antibiotics (penicillins and cephalosporins) is located below the sulfonamide line while the line relating to tetracyclines is located above it. When the reagent from the microwell is re-suspended with a milk sample, both receptors and monoclonal antibodies will bind the corresponding analytes if present during the first 3-minute incubation at 40°C. Afterwards, when the dipstick is dipped into the milk, the liquid starts running vertically on the dipstick and passes through capture zones. When the sample is free of antibiotics, a colour development occurs at the specific capture lines, indicating the absence of the targeted analytes in the milk sample. On the contrary, the presence of antibiotics in the sample will not cause the coloured signal to appear at the specific capture lines.

## IV. Composition of kits

Trisensor Milk Kits contain everything needed to perform 96 measurements.

- 12 pots each with 1 strip of 8 reagent microwells and 8 dipsticks;
- 1 Micropipette of 200  $\mu$ l;
- « Positive Standard » containing powder to reconstitute "positive standard" raw milk of PenG, Sulfadimethoxine and Oxytetracycline (See point X for standards reconstitution). For ease of recognition, the Positive Std is stained by a dye to give it a very pale red colour;
- « Negative Standard » containing powder to reconstitute "negative standard" raw milk (See point X for standards reconstitution). The Negative Std is stained by a dye to give it a very pale green colour;
- One instruction sheet.

## V. Additional material needed

- 1 Heatsensor (40°C incubation, refer to the used Heatsensor manual).
- 1 Readsensensor (optional, refer to the Readsensensor manual).

## VI. General Remarks

1. If instrumental reading is chosen for result interpretation, the Readsensensor must be switched on before the analysis (see Readsensensor manual);
2. At reception, store the kit in a dry place and at a low temperature between 2°C and 8°C. Before opening, let the plastic pots reach room temperature and avoid exposure of the product to moisture and light;
3. Avoid using clotted milk with Trisensor;
4. The best temperature to perform the test is 40°C  $\pm$  3°C. Use the "Heatsensor" (or alternatively a water bath). Any other type of incubator is not appropriate to perform the Trisensor assay. (Refer to the used Heatsensor manual for setting at right temperature and timing);
5. After the second incubation, read the result within a 10-minutes time-frame. Do not attempt to interpret the result after 10 minutes;
6. When drying, the colour intensities of the lines will become sharper;
7. When a positive result is recorded, the test result should be confirmed;
8. Empty one pot before opening another pot.

## VII. Milk powder dilution

In an appropriate bottle, mix 10 g of milk powder with 90 ml of warm (40°C) and distilled water. For an optimal dilution, shake vigorously.

## VIII. Directions for use

*This procedure is described to easily run one single sample or a set of many samples. In that case, try to perform the test in cascade and avoid any delays when mixing reagent and milk but also when adding and removing dipsticks. Make sure you have the same incubation time and temperature for each sample. You shouldn't test more than 8 samples at one time and, if there are more than 3 samples, you should use a multipipette. With more than 8 samples we recommend to share series of maximum 8 samples.*

1. Choose a clean and dry place to perform the test and wash and dry your hands before starting;
2. Connect the Heatsensor (refer to the used Heatsensor manual) and wait until the temperature has stabilised at 40°C.
  - It takes about 10 minutes for the temperature to stabilize at 40°C;
3. Before opening the reagents, take the kit out of the fridge and wait until the temperature of the reagents reaches the ambient temperature. Meanwhile, read the directions for use attentively;
  - There are two main components used for the test which are: dipsticks and the freeze-dried reagents in microwells. Both are stored in the white plastic pots.
4. Determine how many samples are to be tested and write on each tube an identification number;
  - The milk sample must be liquid and homogeneous. There can be neither clots nor sedimentation phases. The ideal temperature of the milk sample is between 4 and 20°C.
5. Open one plastic pot and take out as many microwells as there are milk sample to be tested (the standard positive and negative milk included if necessary);
  - To open a pot of dipsticks, take off the safety ring by pressing it down the pot, take off the ring and get the stopper off the pot with your thumb;
  - The pot with dipsticks should always be well closed after reagents have been taken out;
  - A pot with dipsticks should be emptied before another is opened;
  - Be careful, if you do not intend to use all the 8 microwells, leave the set of 8 caps on the unused ones and do not tear off the strip of the eight caps but leave it on the microwells that will not be used. Do not try to separate individual caps and put them immediately back into the white pot without damaging the dipsticks, close and make sure it is tightly sealed.
6. Place the microwell(s) in the heating block which shows 40°C;
7. Place a new tip on the micropipette and transfer 200 µl of milk into each of the microwells;
8. **Warning: when reagents and milk are in contact, the reaction begins. Mix quickly AND IMMEDIATELY push the START(RUN) button. The 3-minutes countdown starts;**
  - During that first incubation, both receptors and monoclonal antibodies detects whether or not there are antibiotics in the milk sample. It takes 3 minutes for the reaction to be completed.
9. During that time, open the same pot as before, take out as many dipsticks as there are analyses in progress and close the pot. Lay the dipsticks on a clean sheet and write down the number that matches the one of the milk sample;
10. When the 3 minutes are over, i.e. after the sound-signal, press START (STOP)\* again to stop the ringing tone and dip the corresponding dipstick into each of the microwells laid in the incubator;
  - Make sure you put into the corresponding microwell the dipstick with the same identification number as the one of the concerned milk sample;
11. Start the timer for the second incubation by pressing the START (RUN)\* button. The 3-minutes countdown starts;
12. When the 3 minutes are over, i.e. after the sound-signal, press START (STOP)\* again to stop the ringing tone and take the dipsticks out of the microwells to lay them down on a sheet of paper;
13. If you are not planning to perform any other test within the day with Trisensor, put everything back into the box and store it in a fridge at a temperature ranging from 2 to 8°C.  
*\*(Refer to the used Heatsensor manual – For DUO Heatsensor users, steps 9 to 12 are simplified)*

## IX. Visual interpretation of results (see Figure A)

Make visual readings as follows:

14. First check whether the top control line is present. If it is not, regard the analysis as invalid and do not start (or continue) any interpretation.
15. When the top control line can be seen, interpret the three test lines as follows:
16. Examine one test line at a time and compare the intensity of the line colour of the test line with the intensity of the line colour of the control line. (Start with the bottom line of  $\beta$ -lactams for example);
  - If the test line is darker in colour than the control line, the result is **NEGATIVE**, which means that, given the sensitivity of our test, the milk sample contains no antibiotics or antibiotics at a lower level than the value stated in the enclosed table A;

- If the test line is as distinct as or lighter in colour than the control line, the result is POSITIVE, which means that, given the sensitivity of our test, the milk sample contains antibiotics at or above the detection value stated in the enclosed table A.
17. When you have interpreted one test line, do the same for the other lines;
  18. If you hesitate, regard the sample as POSITIVE and confirm your interpretation by performing a second visual reading within the next 15 minutes;
  19. Write down your assessment on each of the dipsticks;
  20. Dipsticks can be archived as a permanent record if required by removing the sample pad and allowing to dry before storage. N.B. Line colour intensity will darken on drying.

If you have a Readsensur, you should read the dipstick within 10 minutes of performing the test. See the user manual of Readsensur.

## X. Negative and Positive standards reconstitution

- Remove the cap and rubber stopper from the vial and add 1 ml of distilled or deionised water;
- Replace the stopper and leave to stand for 1 minute and a half;
- Vigorously shake the sample during 10 seconds;
- Check that the sample is properly dissolved.

The test material is now ready to use.

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